


Anti-Cyclin A2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20248

Item Number: CM0026627



Introduction

Rabbit polyclonal antibody targeting Cyclin A2 (CCNA2), validated for Western blot, IHC-P, IF/ICC, and ELISA in human, mouse, rat samples. Ideal for cell cycle and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	Cyclin A2 Rabbit pAb
Remarks/Alias	CCN1; CCNA; Cyclin A2
Species	Human
GeneID (Human)	890
GeneID	890
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 350-432 of human Cyclin A2 (NP_001228.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20248
Protein Weight	49kDa
Shipping	Ice bag