


Anti-Cyclin A2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P20248

Item Number: CM0003881



Introduction

Rabbit monoclonal antibody targeting human Cyclin A2 (P20248). Validated for WB, IHC-P, IF/ICC, IP, and ELISA; cross-reacts with mouse and rat. Ideal for cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin A2 Rabbit mAb
Remarks/Alias	CCN1; CCNA; Cyclin A2
Species	Human
GeneID (human)	890
GeneID	890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-210 of human Cyclin A2 (P20248)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20248
Protein Weight	49kDa
Shipping	Ice bag