


Anti-Cyclin B1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P14635

Item Number: CM0021183



Introduction

Rabbit polyclonal antibody against human Cyclin B1, validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat CCNB1, a key cell cycle regulator. Suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	Cyclin B1 Rabbit pAb
Remarks/Alias	CCNB; Cyclin B1
Species	Human
GeneID (Human)	891
GeneID	891
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 334-433 of human Cyclin B1 (NP_114172.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14635
Protein Weight	48kDa
Shipping	Ice bag