

Anti-Cyclin B2 Polyclonal Antibody For Wb, Elisa - 095067

Item Number: CM0020042

Introduction

Rabbit polyclonal antibody against human Cyclin B2 (G2/mitotic-specific cyclin-B2). Validated for WB and ELISA, reacts with human and mouse. Ideal for cell cycle and mitosis research.

[Learn More](#)

Parameter	Description
Product Name	Cyclin B2 Rabbit pAb
Remarks/Alias	HsT17299; Cyclin B2
Species	Human
Gene ID (Human)	9133
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Cyclin B2 (NP_004692.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95067
Protein Weight	45kDa
Shipping	Ice bag