

## Anti-Cyclin B2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95067

Item Number: CM0011689

### Introduction

Cyclin B2 Rabbit Monoclonal Antibody (UniProt: O95067) validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Ideal for cell cycle and mitosis research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cyclin B2 Rabbit mAb                                                                                                    |
| Remarks/Alias    | HsT17299; Cyclin B2                                                                                                     |
| Species          | Human                                                                                                                   |
| GeneID (human)   | 9133                                                                                                                    |
| GeneID           | 9133                                                                                                                    |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cyclin B2 (O95067). |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | O95067                                                                                                                  |
| Protein Weight   | 45kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |