

Anti-Cyclin B2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95067

Item Number: CM0029486

Introduction

Cyclin B2 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA applications. Detects human, mouse, and rat Cyclin B2, a key cell cycle regulator. Ideal for mitosis and cell division research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin B2 Rabbit pAb
Remarks/Alias	HsT17299; Cyclin B2
Species	Human
GenelD (human)	9133
GenelD	9133
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95067
Protein Weight	45kDa
Shipping	Ice bag