

Anti-Cyclin C Polyclonal Antibody For Wb And Elisa - P24863

Item Number: CM0015046



Introduction

Rabbit polyclonal antibody against human Cyclin C, ideal for Western blot and ELISA applications. Cross-reacts with mouse and rat to support transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin C Rabbit pAb
Remarks/Alias	CycC; SRB11; hSRB11; Cyclin C
Species	Human
GeneID (Human)	892
GeneID	892
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human Cyclin C (NP_005181.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24863
Protein Weight	33kDa
Shipping	Ice bag