

Anti-Cyclin C Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24863

Item Number: CM0030763



Introduction

Rabbit polyclonal antibody targeting human Cyclin C (CCNC, hSRB11), validated for WB, IHC-P, and ELISA. Detects endogenous levels of Cyclin C in human, mouse, and rat. Applications include transcription and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin C Rabbit pAb
Remarks/Alias	CycC, SRB11, hSRB11, Cyclin C
Species	Human
GeneID (human)	892
GeneID	892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human Cyclin C (NP_005181.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24863
Protein Weight	33kDa
Shipping	Ice bag