


Anti-Cyclin D1 Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P24385

Item Number: CM0016698



Introduction

High-quality rabbit polyclonal antibody against human Cyclin D1 (BCL1). Validated for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin D1 Rabbit pAb
Remarks/Alias	BCL1; PRAD1; U21B31; D11S287E; Cyclin D1
Species	Human
GeneID (Human)	595
GeneID	595
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-295 of human Cyclin D1 (NP_444284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24385
Protein Weight	34kDa
Shipping	Ice bag