


Anti-Cyclin D1 Rabbit Polyclonal Antibody For Wb - P24385

Item Number: CM0016770



Introduction

High-quality rabbit polyclonal antibody targeting Cyclin D1 (CCND1, BCL1, PRAD1). Validated for WB, IF/ICC, IP, and ELISA applications; cross-reacts with human, mouse, and rat. Ideal for investigating G1/S cell cycle regulation and oncoprotein signaling.

[Learn More](#)

Parameter	Value
Product Name	Cyclin D1 Rabbit pAb
Remarks/Alias	BCL1; PRAD1; U21B31; D11S287E; Cyclin D1
Species	Human
GeneID (Human)	595
GeneID	595
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-295 of human Cyclin D1 (NP_444284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24385
Protein Weight	34kDa
Shipping	Ice bag