

Anti-Cyclin D2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30279

Item Number: CM0016783



Introduction

High-quality anti-Cyclin D2 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat Cyclin D2 protein (MW 33 kDa). Suitable for cell cycle and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Cyclin D2 Rabbit pAb
Remarks/Alias	MPPH3; KIAK0002; Cyclin D2
Species	Human
GeneID (Human)	894
GeneID	894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human Cyclin D2 (NP_001750.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30279
Protein Weight	33kDa
Shipping	Ice bag