


Anti-Cyclin D3 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P30281

Item Number: CM0016557



Introduction

Rabbit polyclonal antibody against human G1/S-specific cyclin-D3 (CCND3, SWISS P30281). Validated by knockout, suitable for WB and ELISA. Cross-reacts with mouse for cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cyclin D3 Rabbit pAb
Remarks/Alias	CCND3; cyclin D3; Cyclin D3; D3
Species	Human
GeneID (human)	896
GeneID	896
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30281
Protein Weight	33kDa
Shipping	Ice bag