

Anti-Cyclin E1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P24864

Item Number: CM0005144



Introduction

Rabbit monoclonal antibody against human Cyclin E1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Designed for cell cycle research.

[Learn More](#)

Parameter	Specification
Product Name	Cyclin E1 Rabbit PolymAb®
Remarks/Alias	CCNE; pCCNE1
Species	Human
GeneID (Human)	898
GeneID	898
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-400 of human Cyclin E1 (NP_001229.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24864
Protein Weight	47kDa
Shipping	Ice bag