

Anti-Cyclin E1 Rabbit Monoclonal Antibody For Wb & Elisa - P24864

Item Number: CM0006947



Introduction

Anti-Cyclin E1 rabbit monoclonal antibody validated for WB and ELISA in human, mouse, and rat. Targets the G1/S-specific cyclin-E1, a key cell cycle regulator. Recombinant monoclonal; immunogen: synthetic peptide (aa 300-400). 47kDa, nucleus.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E1 Rabbit mAb
Remarks/Alias	CCNE; pCCNE1; Cyclin E1
Species	Human
Gene ID (Human)	898
Gene ID	898
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human Cyclin E1 (NP_001229.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24864
Protein Weight	47kDa
Shipping	Ice bag