


Anti-Cyclin E2 Monoclonal Antibody For Wb, Elisa - O96020

Item Number: CM0007995



Introduction

Cyclin E2 Rabbit Monoclonal Antibody validated for WB and ELISA. Recognizes human Cyclin E2, a key regulator of G1/S transition, with a molecular weight of ~47 kDa. Ideal for cell cycle, proliferation, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E2 Rabbit mAb
Remarks/Alias	CYCE2; Cyclin E2
Species	Human
Gene ID (Human)	9134
Gene ID	9134
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Cyclin E2 (O96020).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O96020
Protein Weight	47kDa
Shipping	Ice bag