


Anti-Cyclin E2 Polyclonal Antibody For Wb, If/Icc, Elisa - 096020

Item Number: CM0029666



Introduction

Rabbit polyclonal antibody against human Cyclin E2 (CCNE2), suitable for Western blot, immunofluorescence, and ELISA. Detects a ~47 kDa band. Ideal for cell cycle and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E2 Rabbit pAb
Remarks/Alias	CYCE2; Cyclin E2
Species	Human
GeneID (Human)	9134
GeneID	9134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-404 of human Cyclin E2 (NP_477097.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	096020
Protein Weight	47kDa
Shipping	Ice bag