


Anti-Cyclin E2 Polyclonal Antibody For Wb, Ihc-P, Elisa - O96020

Item Number: CM0027439



Introduction

Rabbit polyclonal antibody against human Cyclin E2 (CYCE2). Suitable for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for cell cycle, cancer, and proliferation research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin E2 Rabbit pAb
Remarks/Alias	CYCE2; Cyclin E2
Species	Human
Gene ID (Human)	9134
Gene ID	9134
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 255-404 of human Cyclin E2 (NP_477097.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96020
Protein Weight	47kDa
Shipping	Ice bag