

Anti-Cyclin G1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P51959

Item Number: CM0026178

Introduction

Rabbit polyclonal antibody against Cyclin G1 (CCNG1), validated for WB, IHC-P, ELISA in human, mouse, and rat. Cyclin G1 is associated with G2/M phase arrest in DNA damage response. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin G1 Rabbit pAb
Remarks/Alias	CCNG; Cyclin G1
Species	Human
Gene ID	900 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-295 of human Cyclin G1 (NP_004051.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51959
Protein Weight	34kDa
Shipping	Ice bag