


Anti-Cyclin H Monoclonal Antibody For Wb, Elisa - P51946

Item Number: CM0002992



Introduction

Rabbit monoclonal antibody against Cyclin H (CAK) for WB, ELISA applications. Reacts with human and mouse. Involved in cell cycle regulation and transcription via CDK-activating kinase complex.

[Learn More](#)

Parameter	Value
Product Name	Cyclin H Rabbit mAb
Remarks/Alias	CAK; p34; p37; Cych; Cyclin H
Species	Human
GeneID	902
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cyclin H (P51946)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51946
Protein Weight	38kDa
Shipping	Ice bag