


## Anti-Cyclin H Polyclonal Antibody For Wb, Ip, Elisa - P51946

Item Number: CM0017669



### Introduction

Rabbit polyclonal antibody targeting human Cyclin H, a CAK complex subunit. Validated for WB, IP, and ELISA; supports cell cycle and transcription regulation studies. (38 kDa)

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cyclin H Rabbit pAb                                                                                                  |
| Remarks/Alias    | CAK; p34; p37; Cych; Cyclin H                                                                                        |
| Species          | Human                                                                                                                |
| Gene ID          | 902                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human Cyclin H (NP_001230.1). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IP, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                                |
| SWISS            | P51946                                                                                                               |
| Protein Weight   | 38kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |