


Anti-Cyclin H Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51946

Item Number: CM0011241



Introduction

Cyclin H Rabbit mAb is a rabbit monoclonal antibody recognizing human and mouse Cyclin H, intended for WB, IHC-P, IF/ICC, and ELISA. Validated for research use in cell cycle and transcription regulation studies. 38 kDa.

[Learn More](#)

Parameter	Value
Product Name	Cyclin H Rabbit mAb
Remarks/Alias	CAK; p34; p37; Cych; Cyclin H
Species	Human
Gene ID (Human)	902
Gene ID	902
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 224-323 of human Cyclin H (P51946).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51946
Protein Weight	38kDa
Shipping	Ice bag