


Anti-Cyclin H Rabbit Polyclonal Antibody For Wb, Elisa - P51946

Item Number: CM0020940



Introduction

Cyclin H Rabbit Polyclonal Antibody (P51946) validated for Western blot and ELISA. Targets human Cyclin H, a component of the CDK-activating kinase complex, essential for cell cycle regulation and RNA polymerase II transcription. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Cyclin H Rabbit pAb
Remarks/Alias	CAK; p34; p37; Cych; Cyclin H
Species	Human
GeneID (Human)	902
GeneID	902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human Cyclin H (NP_001230.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51946
Protein Weight	38kDa
Shipping	Ice bag