


Anti-Cyclin T1 Polyclonal Antibody For Wb, Elisa - O60563

Item Number: CM0023261



Introduction

High-quality rabbit polyclonal antibody targeting human Cyclin T1 (CCNT/CYCT1), validated for WB and ELISA. Ideal for studying transcriptional elongation, P-TEFb complex, and HIV Tat interaction.

[Learn More](#)

Parameter	Details
Product Name	Cyclin T1 Rabbit pAb
Remarks/Alias	CCNT; CYCT1; HIVE1; Cyclin T1
Species	Human
GeneID (Human)	904
GeneID	904
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 220-320 of human Cyclin T1 (NP_001231.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60563
Protein Weight	81kDa
Shipping	Ice bag