

Anti-Cyclophilin 40 Monoclonal Antibody For Wb - Q08752

Item Number: CM0009033

Introduction

Rabbit monoclonal antibody specific for Cyclophilin 40/PPID. Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Ideal for protein folding, HSP90 chaperone, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Cyclophilin 40 Rabbit mAb
Remarks/Alias	CYPD; CYP-40; Cyclophilin 40
Species	Human
Gene ID (Human)	5481
Gene ID	5481
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 271-370 of human Cyclophilin 40 (Q08752).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08752
Protein Weight	41kDa
Shipping	Ice bag