

Anti-Cyclophilin 40 Polyclonal Antibody For Wb, If/Icc, Elisa - Q08752

Item Number: CM0027901

Introduction

Rabbit polyclonal antibody against human Cyclophilin 40 (CYP-40), validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat; UniProt Q08752.

[Learn More](#)

Parameter	Value
Product Name	Cyclophilin 40 Rabbit pAb
Remarks/Alias	CYPD; CYP-40; Cyclophilin 40
Species	Human
Gene ID (Human)	5481
Gene ID	5481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 161-370 of human Cyclophilin 40 (NP_005029.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08752
Protein Weight	41kDa
Shipping	Ice bag