

Anti-Cyclophilin B Monoclonal Antibody For Wb, Elisa - P23284

Item Number: CM0012528



Introduction

Rabbit monoclonal antibody targeting human Cyclophilin B (PPIB, CYPB). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for protein folding and chaperone function studies.

[Learn More](#)

Parameter	Value
Product Name	Cyclophilin B Rabbit mAb
Remarks/Alias	OI9; CYPB; SCYLP; CYP-S1; HEL-S-39; Cyclophilin B
Species	Human
Gene ID (Human)	5479
Gene ID	5479
Immunogen	Synthetic peptide corresponding to amino acids 100-200 of human Cyclophilin B (P23284)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23284
Protein Weight	24kDa
Shipping	Ice bag