

Anti-Cyfp1 Polyclonal Antibody For Wb, Ip, Elisa - Q7L576

Item Number: CM0015892



Introduction

Rabbit polyclonal antibody against human CYFIP1, suitable for WB, IP, ELISA. Specific immunogen is N-terminal amino acids 1-270. Cross-reacts with human and mouse. Ideal for studying translational repression, actin remodeling, and neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	CYFIP1 Rabbit pAb
Remarks/Alias	SHYC; SRA1; SRA-1; P140SRA-1; CYFIP1
Species	Human
GeneID (human)	23191
GeneID	23191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CYFIP1 (NP_055423.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7L576
Protein Weight	145kDa
Shipping	Ice bag