

Anti-Cyfp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q7L576

Item Number: CM0003146



Introduction

CYFIP1 Rabbit mAb (Q7L576) validated for WB, IF/ICC, ELISA; reacts with Human, Mouse, Rat. Suitable for studies on translational repression, actin remodeling, and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	CYFIP1 Rabbit mAb
Remarks/Alias	SHYC; SRA1; SRA-1; P140SRA-1; CYFIP1
Species	Human
GeneID	23191
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYFIP1 (Q7L576)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L576
Protein Weight	145kDa
Shipping	Ice bag