

Anti-Cyld Polyclonal Antibody For If/Icc, If-P, Elisa - Q9Nqc7

Item Number: CM0022874



Introduction

Rabbit polyclonal antibody against human CYLD (Ubiquitin carboxyl-terminal hydrolase CYLD). Validated for IF/ICC, IF-P, ELISA. Reacts with human and mouse. Suitable for research on NF- κ B and Wnt signaling, microtubule dynamics, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	CYLD Rabbit pAb
Remarks/Alias	EAC; MFT; SBS; TEM; BRSS; CDMT; MFT1; CYLD1; CYLDI; USPL2; FTDALS8
Species	Human
Gene ID (Human)	1540
Gene ID	1540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human CYLD (NP_056062.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQC7
Protein Weight	107kDa
Shipping	Ice bag