

Anti-Cyld Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nqc7

Item Number: CM0030489



Introduction

CamelBio's CYLD Rabbit pAb (Q9NQC7) is a high-quality polyclonal antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for NF-κB and Wnt signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CYLD Rabbit pAb
Remarks/Alias	EAC; MFT; SBS; TEM; BRSS; CDMT; MFT1; CYLD1; CYLDI; USPL2; FTDALS8; CYLD
Species	Human
GeneID (human)	1540
GeneID	1540
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQC7
Protein Weight	107kDa
Shipping	Ice bag