

Anti-Cynomolgus Monkey IgD Rabbit Monoclonal Antibody - Monkey IgD

Item Number: CM0008682



Introduction

Rabbit monoclonal antibody specific for cynomolgus monkey immunoglobulin delta heavy chain. Validated for IHC-P, FC, ELISA. For immune profiling and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Monkey IgD Rabbit mAb
Species	Cynomolgus
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC, ELISA
Cross Reactivity	Cynomolgus
Shipping	Ice bag