

## Anti-Cyp11A1 Rabbit Monoclonal Antibody For Wb And Elisa - P05108

Item Number: CM0004696



### Introduction

Rabbit monoclonal antibody against human CYP11A1, validated for WB and ELISA on human and mouse samples. Targets mitochondrial cholesterol side-chain cleavage enzyme involved in steroidogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | CYP11A1 Rabbit mAb                                                                                                  |
| Remarks/Alias    | CYP11A1; CYP11A; CYPXIA1; P450SCC; cytochrome P450 family 11 subfamily A member 1                                   |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 1583                                                                                                                |
| Gene ID          | 1583                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 40-320 of human CYP11A1 (NP_000772.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | P05108                                                                                                              |
| Protein Weight   | 42kDa/60kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                             |