

Anti-Cyp11A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P05108

Item Number: CM0019452



Introduction

CYP11A1 rabbit polyclonal antibody validated for Western blot, immunofluorescence, IHC, and ELISA. Detects human, mouse, rat cholesterol side-chain cleavage enzyme for steroidogenesis research. - P05108

[Learn More](#)

Parameter	Value
Product Name	CYP11A1 Rabbit pAb
Remarks/Alias	CYP11A; CYPXIA1; P450SCC; CYP11A1
Species	Human
GeneID (Human)	1583
GeneID	1583
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 40-320 of human CYP11A1 (NP_000772.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05108
Protein Weight	60kDa
Shipping	Ice bag