


Anti-Cyp11B1 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P15538

Item Number: CM0017366



Introduction

Rabbit polyclonal antibody targeting human CYP11B1, validated for IHC-P, IF/ICC, and ELISA. Ideal for investigating glucocorticoid biosynthesis, adrenal cortex function, and steroidogenic pathways. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CYP11B1 Rabbit pAb
Remarks/Alias	FHI; CPN1; CYP11B; P450C11; CYP11B1
Species	Human
GeneID (Human)	1584
GeneID	1584
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human CYP11B1 (NP_000488.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P15538
Protein Weight	58kDa
Shipping	Ice bag