

Anti-Cyp17A1 Rabbit Polyclonal Antibody For Wb, Elisa - P05093

Item Number: CM0017329



Introduction

CYP17A1 rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human CYP17A1 (57kDa), cross-reacts with rat. Ideal for studying steroid hormone biosynthesis and cytochrome P450 function. UniProt P05093.

[Learn More](#)

Parameter	Value
Product Name	CYP17A1 Rabbit pAb
Remarks/Alias	CPT7; CYP17; S17AH; P450C17; CYP17A1
Species	Human
Gene ID (Human)	1586
Gene ID	1586
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-508 of human CYP17A1 (NP_000093.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P05093
Protein Weight	57kDa
Shipping	Ice bag