

Anti-Cyp1A2 Polyclonal Antibody For Wb, If-P, Elisa - P05177

Item Number: CM0014740



Introduction

Rabbit polyclonal CYP1A2 antibody (SWISS P05177) validated for WB, IF-P, ELISA; cross-reacts with human, mouse, rat. Ideal for studying cytochrome P450 drug metabolism, lipid signaling, and toxicology pathways.

[Learn More](#)

Parameter	Value
Product Name	CYP1A2 Rabbit pAb
Remarks/Alias	CP12; CYP1A2; P3-450; P450(PA); CYP1A2
Species	Human
GeneID (Human)	1544
GeneID	1544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 205-305 of human CYP1A2 (NP_000752.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05177
Protein Weight	58kDa
Shipping	Ice bag