


Anti-Cyp1B1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q16678

Item Number: CM0018501



Introduction

Rabbit anti-human CYP1B1 polyclonal antibody validated for Western Blot, IF/ICC, and ELISA. Targets cytochrome P450 1B1 involved in steroid and xenobiotic metabolism, suitable for disease research and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	CYP1B1 Rabbit pAb
Remarks/Alias	CP1B; ASGD6; GLC3A; CYP1B1; P4501B1; CYP1B1
Species	Human
GeneID (Human)	1545
GeneID	1545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-543 of human CYP1B1 (NP_000095.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q16678
Protein Weight	61kDa
Shipping	Ice bag