


Anti-Cyp1B1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16678

Item Number: CM0014164



Introduction

This rabbit polyclonal antibody against human CYP1B1 is suitable for Western blot, immunofluorescence/immunocytochemistry, and ELISA. It cross-reacts with mouse and rat samples for cytochrome P450 metabolism and disease research.

[Learn More](#)

Parameter	Value
Product Name	CYP1B1 Rabbit pAb
Remarks/Alias	CP1B; ASGD6; GLC3A; CYP1B1; P4501B1; CYP1B1
Species	Human
GeneID (Human)	1545
GeneID	1545
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CYP1B1 (NP_000095.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16678
Protein Weight	61kDa
Shipping	Ice bag