

Anti-Cyp27A1 Polyclonal Antibody For Wb, Elisa - Q02318

Item Number: CM0021260

Introduction

CYP27A1 Rabbit pAb validated for WB, ELISA; reacts with human, mouse, rat. Ideal for cholesterol metabolism, bile acid synthesis, and related disease research.

[Learn More](#)

Parameter	Specification
Product Name	CYP27A1 Rabbit pAb
Remarks/Alias	CTX; CP27; CYP27; CYP27A1
Species	Human
Gene ID (Human)	1593
Gene ID	1593
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human CYP27A1 (NP_000775.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02318
Protein Weight	60kDa
Shipping	Ice bag