

Anti-Cyp27A1 Rabbit Monoclonal Antibody For Wb And Elisa - Q02318

Item Number: CM0007263

Introduction

Recombinant rabbit monoclonal antibody targeting CYP27A1, sterol 26-hydroxylase. Validated for WB and ELISA, cross-reacts with human, mouse, rat. For studying cholesterol metabolism and bile acid synthesis. Ideal for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CYP27A1 Rabbit mAb
Remarks/Alias	CTX; CP27; CYP27; CYP27A1
Species	Human
GeneID (human)	1593
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 431-530 of human CYP27A1 (NP_000775.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q02318
Protein Weight	60kDa
Shipping	Ice bag