


Anti-Cyp2C18 Rabbit Polyclonal Antibody For Wb, Elisa - P33260

Item Number: CM0029936



Introduction

Rabbit polyclonal antibody against human CYP2C18, validated for WB and ELISA, cross-reactive with mouse and rat. Suitable for retinoid metabolism and cytochrome P450 studies. Immunogen: amino acids 231-490.

[Learn More](#)

Parameter	Value
Product Name	CYP2C18 Rabbit pAb
Remarks/Alias	CPCI; CYP2C; CYP2C17; P450IIC17; P450-6B/29C; CYP2C18
Species	Human
Gene ID (Human)	1562
Gene ID	1562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-490 of human CYP2C18 (NP_000763.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33260
Protein Weight	56kDa
Shipping	Ice bag