


Anti-Cyp2D6 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P10635

Item Number: CM0010564



Introduction

High-specificity rabbit monoclonal antibody against CYP2D6, validated for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Detects ~56 kDa protein. Ideal for drug metabolism and pharmacogenomics studies.

[Learn More](#)

Parameter	Value
Product Name	CYP2D6 Rabbit mAb
Remarks/Alias	CPD6; CYP2D; CYP2DL1; CYP2D6; P450C2D; P450DB1; CYP2D7AP; CYP2D7BP; CYP2D7P2; CYP2D8P2; P450-DB1; CYP2D6
Species	Human
GeneID (human)	1565
GeneID	1565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-497 of human CYP2D6 (P10635)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10635
Protein Weight	56kDa
Shipping	Ice bag