


Anti-Cyp2D6 Rabbit Polyclonal Antibody For Wb And Elisa - P10635

Item Number: CM0014655



Introduction

CYP2D6 Rabbit Polyclonal Antibody validated for WB and ELISA; detects human CYP2D6 with cross-reactivity to mouse and rat. Ideal for studying cytochrome P450 drug metabolism and lipid signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	CYP2D6 Rabbit pAb
Remarks/Alias	CPD6; CYP2D; CYP2DL1; CYP2D6; P450C2D; P450DB1; CYP2D7AP; CYP2D7BP; CYP2D7P2; CYP2D8P2; P450-DB1; CYP2D6
Species	Human
GeneID (Human)	1565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of human CYP2D6 (NP_000097.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10635
Protein Weight	56kDa
Shipping	Ice bag