

Anti-Cyp2E1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05181

Item Number: CM0021850



Introduction

Rabbit anti-CYP2E1 polyclonal antibody for research use in Western blot, IF, IHC, and ELISA. Detects human, mouse, and rat CYP2E1, a liver cytochrome P450 involved in fatty acid and xenobiotic metabolism.

[Learn More](#)

Parameter	Value
Product Name	CYP2E1 Rabbit pAb
Remarks/Alias	CPE1; CYP2E; P450-J; P450C2E; CYP2E1
Species	Human
GeneID (Human)	1571
GeneID	1571
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 234-493 of human CYP2E1 (NP_000764.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05181
Protein Weight	57 kDa
Shipping	Ice bag