

Anti-Cyp2F1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P24903

Item Number: CM0030024

Introduction

Rabbit polyclonal antibody against CYP2F1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. CYP2F1 is a cytochrome P450 enzyme involved in pneumotoxicant metabolism; ideal for lung toxicology studies.

[Learn More](#)

Parameter	Value
Product Name	CYP2F1 Rabbit pAb
Remarks/Alias	C2F1; CYP2F; CYP11F1; CYP2F1
Species	Human
Gene ID (Human)	1572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-290 of human CYP2F1 (NP_000765.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24903
Protein Weight	56kDa
Shipping	Ice bag