

Anti-Cyp2R1 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Vvx0

Item Number: CM0014329



Introduction

Rabbit polyclonal antibody to human CYP2R1 (Vitamin D 25-hydroxylase), validated for Western blot and ELISA. Cross-reacts with mouse and rat; detects ~57 kDa protein. Ideal for vitamin D metabolism research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CYP2R1 Rabbit pAb
Remarks/Alias	CYP2R1
Species	Human
GeneID (Human)	120227
GeneID	120227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 252-501 of human CYP2R1 (NP_078790.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6VVX0
Protein Weight	57kDa
Shipping	Ice bag