

Anti-Cyp39A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nyl5

Item Number: CM0022964

Introduction

Rabbit polyclonal antibody targeting human CYP39A1 (24-hydroxycholesterol 7-alpha-hydroxylase). Validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for bile acid synthesis and cholesterol metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CYP39A1 Rabbit pAb
Remarks/Alias	CYP39A1
Species	Human
GeneID (Human)	51302
GeneID	51302
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYP39A1 (NP_057677.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYL5
Protein Weight	54kDa
Shipping	Ice bag