

Anti-Cyp39A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nyl5

Item Number: CM0021331

Introduction

High-quality rabbit polyclonal antibody against human CYP39A1, validated for WB and ELISA. Detects mouse and rat orthologs. Suitable for studying neural cholesterol metabolism and liver-specific oxysterol 7-alpha-hydroxylase function.

[Learn More](#)

Parameter	Specification
Product Name	CYP39A1 Rabbit pAb
Species	Human
Gene ID (Human)	51302
Gene ID	51302
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NYL5
Protein Weight	54kDa
Shipping	Ice bag