


Anti-Cyp3A4 Monoclonal Antibody For Wb, If-P, Elisa - P08684

Item Number: CM0006812



Introduction

Rabbit monoclonal antibody targeting human CYP3A4 for Western blot, immunofluorescence, and ELISA applications. Detects a ~57 kDa protein. Ideal for studying drug metabolism and CYP3A4-related pathways.

[Learn More](#)

Parameter	Value
Product Name	CYP3A4 Rabbit mAb
Remarks/Alias	HLP; CP33; CP34; CYP3A; NF-25; VDDR3; CYP3A3; CYP11A3; CYP11A4; P450C3; CYP11A4; P450PCN1; CYP3A4
Species	Human
Gene ID (Human)	1576
Gene ID	1576
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-503 of human CYP3A4 (NP_059488.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P08684
Protein Weight	57kDa
Shipping	Ice bag