

Anti-Cyp3A4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08684

Item Number: CM0026270



Introduction

High-quality rabbit polyclonal anti-CYP3A4 antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Target for drug metabolism, steroid biosynthesis, and vitamin D catabolism research.

[Learn More](#)

Parameter	Value
Product Name	CYP3A4 Rabbit pAb
Remarks/Alias	HLP; CP33; CP34; CYP3A; NF-25; VDDR3; CYP3A3; P450C3; CYP11A3; CYP11A4; P450PCN1; CYP3A4
Species	Human
GeneID (human)	1576
GeneID	1576
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08684
Protein Weight	57kDa
Shipping	Ice bag